




## RECOMMENDED LENGTH OF NEEDLES

	Age	 BMI		Needle length*
	Children < 12 Very slim adults		Without skin fold straight injection	With skin fold, angled or straight injection
		all	4 mm, 5 mm, 6 mm	6 mm
	Teenagers (12-18 years)	average	4 mm, 5 mm, 6 mm	8 mm
		overweight (BMI > 25)	4 mm, 5 mm, 6 mm	8 mm, 10 mm
	Adults	average	4 mm, 5 mm, 6 mm, 8 mm	8 mm
		overweight (BMI > 25)	4 mm, 5 mm, 6 mm, 8 mm	8 mm, 10 mm, 12 mm**

\* Consult your healthcare professional for assistance in choosing the appropriate needle length and injection technique.

\*\* Also recommended for people injecting large dosage or with limited dexterity/strength.

Droplet® pen needles have been tested for compliance with ISO 11608-2:2012.

Droplet® pen needles work with injection pens available on the US market produced by the following companies:

**ELI LILLI, OWEN MUMFORD, HASELMEIER**

**A full list of injection pens working with Droplet® pen needles →**