

COMPATIBILITY STATEMENT

HTL-STREFA S.A. is the manufacturer of personal lancets, type 560, which are marketed in the following versions: 28G, 30G, 33G

This statement is based on successfully passed tests conducted with lancing devices listed in Appendix 1.

Personal lancets type 560 were tested to evaluate their compatibility with certain commercially available lancing devices from various manufacturers.

The evaluation included visual assessment of device features, strength testing, and the measurement of parameters that may affect compatibility between lancet and lancing device.

The conducted tests and evaluation of compatibility were used to establish the current list of lancing devices compatible with personal lancets type 560.

A non-exhaustive list of compatible lancing devices can be found in Appendix 1.

The brand names of pen injector devices are trademarks of their respective owners and may vary across different markets.

Appendix 1

List of lancing devices work with personal lancets type 560 manufactured by HTL-STREFA S.A.

Table 1. Summary of test results on compatibility personal lancet type 560 with lancing devices.

Manufacturer / Distributor	Brand name	Product picture	Test result Assessment whether the personal lancet works with type 560 personal lancets (YES/NO)
HTL-STREFA S.A.	droplet® (type 710)		YES
	droplet® (type 700)		YES
Abbott	FreeStyle		YES
Bayer	Microlet® 2		YES
ADVOCATE	ADVOCATE®		YES
MHC Medical	Easy Touch®		YES
Owen Mumford	Autolet® Plus		YES
Cardinal Health Essentials	Mini Lancing Device		YES
ReliOn	ReliOn™ Lancing Device		YES
Sanofi	BG Star®		YES

Manufacturer / Distributor	Brand name	Product picture	Test result Assessment whether the personal lancet works with type 560 personal lancets (YES/NO)
Simple Diagnostics	Simple Diagnostics®		YES
Ypsomed	my life Softlance		YES
ADC	Microlet® next		YES
Genteel	Genteel® lancing device		YES
PIKDARE S.p.A	Digitest		YES
STERILANCE MEDICAL (SUZHOU) INC.	Digitest		YES
ACON Laboratories, Inc.	On Call® Stealth® Lancing Device		YES
Future Diagnostics LLC	CareTouch™		YES
Huahong	Blood Lancet Pen Painless Sterile Gold Plated		YES
Bionime Corporation GE	Bionime Lancing Device		YES
iHealth	iHealth®		YES
Owen Mumford	Autolet® Impression		YES

Manufacturer / Distributor	Brand name	Product picture	Test result Assessment whether the personal lancet works with type 710 lancing device (YES/NO)
A. Menarini	Glucoject® Dual		YES
	Glucoject® Dual Plus		YES
	MenaLancet® DUAL		YES
Lifescan Johnson & Johnson	Penlet™ Plus		YES
	ONETOUCH		YES
	OneTouch UltraSoft		YES
BD	BD Optimus™		YES
Trividia	TRUEdraw™		YES
ReliaMed®	Mini lancing device		YES
	Pen-style Mini Lancing Device		YES
Rite Aid Pharmacy	Rite Aid®		YES
SteriLance	SteriLance™		YES