

# COMPATIBILITY STATEMENT

HTL-STREFA S.A.  
is the manufacturer  
of lancing device,  
type 710

This statement is based on successfully passed tests conducted with personal lancets listed in Appendix 1.

Lancing device type 710 was tested to evaluate its compatibility with certain commercially available personal lancets from various manufacturers.

The tests consisted of visual assessment of the device features, strength measurements, and dimensional measurements that could affect the device's compatibility with personal lancets.

The conducted tests and compatibility evaluation were used to confirm the list of personal lancets compatible with the lancing device type 710 provided below.

### Appendix 1

List of personal lancets working with lancing device type 710 manufactured by HTL-STREFA S.A.

Table 1. Summary of test results for the compatibility of lancing device type 710 with personal lancets.

Manufacturer / Distributor	Brand name	Product picture	Test result Assessment whether the personal lancet works with type 710 lancing device (YES/NO)
HTL-STREFA S.A.	droplet® personal lancet (type 560)		YES
HTL-STREFA S.A.	droplet® blood lancet (type 561)		YES