

# STATEMENT

**HTL-STREFA S.A. is the manufacturer of pen needles type 810 which are marketed in the following versions:**

1. 0.23 mm x 4 mm (32G x 5/32")
2. 0.23 mm x 5 mm (32G x 3/16")
3. 0.23 mm x 6 mm (32G x 1/4")
4. 0.23 mm x 8 mm (32G x 5/16")
5. 0.25 mm x 5 mm (31G x 3/16")
6. 0.25 mm x 6 mm (31G x 1/4")
7. 0.25 mm x 8 mm (31G x 5/16")
8. 0.33 mm x 10 mm (29G x 3/8")
9. 0.33 mm x 12 mm (29G x 1/2")
10. 0.30 mm x 8 mm (30G x 5/16")
11. 0.20 mm x 4 mm (33G x 5/32")
12. 0.18 mm x 3.5 mm (34G x 3.5 mm)

Product availability per brand name and per version (gage & length) may vary on various markets.

Pen needles have been tested for compliance with standard **EN ISO 11608-2:2012 Needle-based injection systems for medical use -- Requirements and test methods -- Part 2: Needles.**

Based on testing in accordance with Clause 4.9 and 11 of EN ISO 11608-2:2012, compatibility with needle-based injection systems (NIS), hereinafter called pen injector devices can be claimed. Pen needles have passed testing required by ISO 11608-2.

A non-exhaustive list of compatible pen injector devices is to be found in the Appendix 1 to this statement (List of pen injector devices compatible with pen needles manufactured by HTL-STREFA).






















Brand names of pen injector devices are trademarks of their respective owners and may vary on different markets.



**Enclosed:**

Appendix 1 - List of pen injector devices compatible with pen needles manufactured by HTL-STREFA

**Table 1 - Summary of test results on compatibility pen needles with needle-based injection system**

Standard		EN ISO 11608-2:2012 clause 4.9 & 11	
Needle-based injection system			Test result
Manufacturer / Distributor	Brand name	Product picture	Compatible / Non-compatible
B. Braun	Omnican® Pen 31		✓
	Omnican® Pen 32		✓
Owen Mumford / Berlin-Chemie	BerliPen® 301		✓
	BerliPen® 302		✓
Haselmeier / Berlin-Chemie	BerliPen® JUNIOR		✓
	BerliPen® areo2		✓
	BerliPen® PRECISION		✓
	BerliPen® areo		✓
Haselmeier / Sanofi-Aventis	JuniorSTAR®		✓
	TactiPen®		✓
Haselmeier	i-pen		✓
	i-pen <sup>®</sup>		✓
Copernicus Sp. z o.o. / Bioton	GensuPen®		✓
Nemera Sp. z o.o. / Bioton	GensuPen 2		✓
	GensuPen IMPROVE		✓
Eli Lilly	Humalog® KwikPen 100 U/ml, 200 U/ml		✓

	HumaPen® Ergo		✓
	HumaPen® Ergo II		✓
	HumaPen® L U X U R A		✓
	HumaPen® L U X U R A HD		✓
	HumaPen® MEMOIR		✓
	HumaPen® savvio		✓
	Humulin® N Pen 100U/ml		✓
	βyetta® Pen		✓
	ABASAGLAR® KwikPen® 100 Units/mL		✓
	Saline KwikPen™		✓
	Forteo® 20 µg/80µl		✓
Novo Nordisk	FlexPen®		✓
	FlexTouch®		✓
	NovoLet®		✓
	NovoPen® 3		✓
	NovoPen® 4		✓
	NovoPen® 5		✓
	NovoPen Echo®		✓
	NovoPen® Junior		✓
	Xultophy®		✓
	Victoza® 6mg/ml		✓

	Tresiba® FlexTouch® 200 units/mL		✓
	Levemir® InnoLet® 100 units/ml		✓
<b>Owen Mumford</b>	Autopen® 24 3.0 ml x 1 unit 1-21 units		✓
	Autopen® 24 3.0 ml x 2 units 2-42 units		✓
	Autopen® Classic 3.0 ml x 1 unit 1-21 units		✓
	Autopen® Classic 3.0 ml x 2 units 2-42units		✓
<b>Sanofi-Aventis</b>	ClikSTAR®		✓
	OptiPen® Pro1		✓
	OptiPen® Pro2		✓
	SoloSTAR®		✓
	Lantus® SoloStar® 100 units/mL		✓
	Lyxumia® 20 µg		✓
<b>Ypsomed</b>	ypsopen®		✓
<b>Ypsomed / Фармстандарт</b>	БиоматикПен®		✓
<b>Ypsomed / Tonghua Dongbao</b>	Gansulin Pen		✓
<b>Ypsomed / Polfa Tarchomin</b>	Polhumin Pen		✓
<b>BD Medical</b>	PuregonPen®		✓
<b>Gani Technology</b>	Xiu Lin Pen		✓
<b>PENDIQ</b>	pendiq		✓
<b>Diamesco Co., Ltd. / pendiq GmbH</b>	pendiq 2.0		✓

**HTL-STREFA S.A.**  
Adamówek 7  
PL-95035 Ozorków, Poland  
T +48 42 270 00 10  
F +48 42 270 00 20  
[www.htl-strefa.pl](http://www.htl-strefa.pl)  
[info@htl-strefa.pl](mailto:info@htl-strefa.pl)

VAT: PL7321880362  
REGON: 472350579  
District Court for Lodz-Strodmiemie in Lodz  
XX Commercial Department of the  
National Court Register  
KRS: 0000256309, share capital:  
2,934,421.10 PLN (paid in full)

5 of 5 ————— 2022-01-24 RA 323a/2019 rev. 3



**NOTE!** The list above is not exhaustive – the brand names may vary on different markets, change due mergers of acquisitions of their owners and private label arrangements.